

LUCIFERASE ASSAY:

1. After cell treatment, wash cells once with PBS. Add **lysis buffer**
200 ul for 6-well plates, 125 ul for 12 well, or 75 ul per well in 24 well plates.
2. Scrape cells free from the well with a rubber policeman or tip of pipette and transfer to a microfuge tube and spin at highest speed for 5 min at RT.
3. Add 50 ul of the cell extract supernatant to 350 ul of freshly-prepared **assay buffer** in an assay tube (12 X 75 mm polypropylene).
4. Place the assay tube in a luminometer and inject with 100 ul of 0.5-1 mM luciferin.
5. Measure peak-integrated luminescence after a 5 sec delay over a 10 sec window.

LYSIS BUFFER

25 mM Tris/phosphate
4 mM EGTA
1% Triton X-100
10% glycerol
2 mM dithiothreitol

For 50 ml:

1.25 ml 1 M Trizma base, 39.25 ml QH₂O, 4 ml 50 mM EGTA. Adjust pH to 7.8 with dilute o-phosphoric acid. Add 5 ml glycerol and 0.5 ml of Triton X-100. Stable if stored at room temperature briefly. Just prior to use 2 uul of 1 M dithiothreitol per 1ml to make complete buffer.

ASSAY BUFFER

25 mM Tris/phosphate
20 mM MgSO₄
4 mM EGTA
2 mM ATP
1 mM dithiothreitol

20 ml

0.5 ml Trizma base, 17.1 ml dH₂O, 0.6 ml 1 M MgSO₄, 1.6 ml 50 mM EGTA. Adjust pH to 7.8 with dilute o-phosphoric acid. Add 200 ul of 200 mM ATP and 20 ul 1 M dithiothreitol just before use.

1 mM D-LUCIFERIN

Directions for a 25 mg sample: Dissolve 0.17 gr dithiothreitol in 110 ml of QH₂O. Add 1 ml of this DTT solution to the 25 mg d-luciferin. Add 20 ul of 10N NaOH to dissolve the luciferin. Dilute this into the remaining 109 ml of DTT solution and store as 1.25 ml aliquots at -20C until use.

REAGENTS

d-luciferin Molecular Proges L2911
ATP Sigma A5394
Triton X-100 Sigma T6878
EGTA Sigma E3889