

Vent PCR

Vent polymerase has exonuclease activity. Because of this, it has a tendency to backup and repair mistaken incorporation. Taq polymerase lacks this activity. For this reason, the fidelity of DNA replication by Vent is much higher than that of Taq. As a general rule, try to amplify up your DNA using Vent if you plan to express it as a protein or otherwise do not wish to tolerate any mutations. However, Vent requires some tooling to work properly, and can be problematic for some fragments sized greater than 1 kb. It all depends....

1. For a standard 100 ul reaction

76.5 – X ul of QH₂O

10 ul of 10X Vent polymerase buffer (supplied with enzyme from NEB)

**X ul of 100 mM MgSO₄

5 ul of Primer A (20 uM stock)

5 ul of Primer B (20 uM stock)

2 ul of dNTP mixture (25 mM stock)

1 ul of template DNA (~1-3 ng/ml stock)

0.5 ul of Vent polymerase (2 U/ul stock)

(** use between 0.5 to 3 ul, see note 4 below)

2. Overlay with a drop of mineral oil

3. PCR cycling

95C for 30 seconds

Primer melting temp minus 4 or 5C for 45 seconds

72C (as a general guideline, use ~ 30 seconds per 500 bp of amplified fragment.

Perform about 30 cycles

4. Note. Sometimes, you can get better results, and even results at all by varying the reaction MgSO₄. NEB supplies 100 mM MgSO₄ for this purpose. Use from 0.5 to 3 ul of this stock MgSO₄ per reaction. It is entirely empirical whether and how much will enhance the product.