

**Protocol for the preparation of dNTP samples  
from mammalian cells**

*For the dNTP assay, prepare extracts from  $2 \times 10^6$  cells using a protocol described below based on the cell type.*

**Non-adherent cells**

- 1) Count cells and prepare cell pellets with  $2 \times 10^6$  cells.
- 2) Wash cell pellet with PBS, and remove as much PBS as possible.
- 3) Lysis the cell pellet by quickly resuspending the cells in ice cold 65% methanol (100ul per  $1 \times 10^6$  cells).
- 4) Vigorously vortex samples for 2 min.
- 5) Completely lyse the cells by incubating at 95 °C for 3 min (make sure an epi tube lid lock is securely in place to prevent opening during the 95 °C incubation).
- 6) Chill on ice for 1 minute (to prevent burning your hands).
- 7) Centrifuge the tube for 3 min at 14K RPM and transfer the 65% methanol solution to a new label tube with identifiable numbers. Discard the tube with the pellet.
- 8) Speed vacuum dry the samples.
- 9) Ship the dNTP samples on dry ice to the following address below.

**Adherent cells**

- 1) Wash cell monolayer twice with PBS. (Always have an extra well for cell counting for each of the experimental groups).
- 2) Lysis the cell monolayer by quickly adding ice cold 65% methanol (100 ul per  $1 \times 10^6$  cells). Scrape the cells off the plate, THOROUGHLY. Add another 100 ul of ice cold 65% methanol to wash the well to recover all the biomaterial. Place all the material into one tube for that sample. Have the tubes on ice while processing multiple samples.
- 3) Vigorously vortex samples for 2 min.
- 4) Completely lyse the cells by incubating at 95 °C for 3 min (make sure an epi tube lid lock is securely in place to prevent opening during the 95 °C incubation).

- 5) Chill on ice for 1 minute (to prevent burning your hands).
- 6) Centrifuge the tube for 3 min at 14K RPM and transfer the 65% methanol solution to a new label tube with identifiable numbers. Discard the tube with the pellet.
- 7) Speed vacuum dry the samples.
- 8) Ship the dNTP samples on dry ice to the following address below.

*FedEx package to this address:*

*Dr. Baek Kim  
c/o Caitlin Shepard  
1760 Haygood Drive  
HSRB, Room E450  
Atlanta, GA 30322  
Phone: (404) 727-1454  
e-mail: baek.kim@emory.edu*